


# The GreenMAMBA retraction made easy.

Product Information:

- 1 deployment in 4 turns
- 2 validated disinfection & sterilisation process
- 3 non reflective shaft to prevent flares
- 4 8cm (straight and with 45° angle) and 12cm wide
- 5 5mm diameter, 33cm length



REF	GMR12	12 cm wide
REF	GMR8B	8 cm wide, 45° angle
REF	GMR8	8 cm wide
REF	FCXL	Flush cylinder for cleaning of the GMR12
REF	FCL	Flush cylinder long for cleaning of the GMR8(B)

2020-08 EN