

# OMNIFinger™

Articulating endoscopic  
clip applier



Easy access and excellent visibility - **60°** of total angulation



Inspired by Diego Gonzalez Rivas

PROVIDING SPECIALISTS WITH  
**INNOVATIVE**  
HEALTHCARE SOLUTIONS

WITH THE CENTRAL FOCUS ON  
PATIENTS' SAFETY AND USERS' COMFORT



# OMNIFinger™

## Articulating endoscopic clip applier

Easy access and excellent visibility - **60°** of total angulation



### KEY FEATURES:

- compatible with non-absorbable Click'aV® and 2<sup>nd</sup> generation Click'aV Plus™ ligating clips,
- 50° articulation + 10° jaws' tip angle gives 60° of total angulation,
- articulation and 360° shaft rotation by index finger allows one-hand operation
- jaws activated by front handle for better stability,
- flushing port facilitates easy cleaning, thus eliminating potential contamination,
- non-detachable design for easier maintenance and durability,
- made of high quality stainless steel,
- colour coded size matching cartridge hue,
- by default without lock, version with lock available for special order,
- PATENT PENDING



**Dr. Diego Gonzalez Rivas**

Degree in Medicine and Surgery by the Universidad de Santiago de Compostela (1992-1998). Resident at the Thoracic Surgery Unit at Juan Canalejo hospital in La Coruna from 1999 to 2004. Medical doctor in Thoracic Surgery since July 2004, with medical practice at the Thoracic Surgery Unit at Complejo Hospitalario Universitario in Santiago de Compostela until April 2005 and from then in the Thoracic Surgery Unit at Complejo Hospitalario Universitario de A Coruna.

**DESIGNED TO ELIMINATE CHALLENGING CASES**



Suitable for use with Click'aV Plus™ (Patent No. 10,265,079) and Click'aV® Ligating Clips



60° of total angulation











One-hand operation



Standard handle: without lock



Optional: handle with lock

REF	SIZE	LOCK	COLOUR	SHAFT DIAMETER (MM)	TROCAR DIAMETER (MM)	WORKING LENGTH (CM)	BOX QUANTITY (PCS)
0301-04MEOMN	M	YES		10	10	36	1
0301-04MEOMNX	M	NO		10	10	36	1
0301-04MLEOMN	ML	YES		10	10	36	1
0301-04MLEOMNX	ML	NO		10	10	36	1
0301-04LEOMN	L	YES		10	10	36	1
0301-04LEOMNX	L	NO		10	10	36	1
0301-04XLEOMN	XL	YES		10	10	36	1
0301-04XLEOMNX	XL	NO		10	10	36	1

Distributed by:

hcp-austria, Aigner Straße 35, A-5026 Salzburg  
Tel.: +43 (0)662/62 12 43, Fax: DW -11  
office@hcp-austria.com, www.hcp-austria.com

*healthcarepartner* GmbH

 **Grena Ltd.**  
1000 Great West Road,  
Brentford, London TW8 9HH,  
United Kingdom

**CUSTOMER SERVICE CENTER**  
NOTTINGHAM OFFICE  
☎ +44 115 9704 800  
☎ +44 115 9704 800  
✉ export@grena.co.uk  
**www.grena.co.uk**



**LONDON CHAMBER**  
COMMERCE *of* AND INDUSTRY