

LogiFlex™



Resposable® Laparoscopic Blunt
Dissection & Gastric Band Placement


For over 30 years Surgical Innovations have been at the forefront of Minimally Invasive Surgery.

As surgical techniques continue to evolve, the demand for solutions which align to those advancements is key, whilst reducing the **impact on the environment**, as societally we become more aware of the footprint the provision of healthcare generates.

Leading the evolution in its peer group, our **hybrid** technology is perfectly positioned to meet these challenges and crucially, is accessible to all.

LogiFlex™ is a dual purpose, **Resposable®** device, for blunt dissection and to assist with gastric band placement. It can also be used for the creation of the retro-gastric tunnel.

LogiFlex™ allows for 360° rotation and 90° articulation, and is supplied as a single use instrument shaft. Fully compatible with our Logic reusable ratcheted handle, minimise the financial and **environmental footprint** generated, without any impact on the quality or performance of the device.



Resposable® Laparoscopic Blunt Dissection & Gastric Band Placement

Ordering Information

Product Code	Product Description
120-7200	Blunt Dissection & Gastric Band Positioning Device (455mm)

➤➤➤ Ratchet Handle

Product Code	Product Description
101-44000	Pistol Grip Handle with Ratchet
101-42000	Horizontal Handle with Ratchet




Leading the Evolution

In Sustainable Surgery.



 Clayton Wood House, 6 Clayton Wood Bank, Leeds, LS16 6QZ, UK

 +44 (0)113 230 7597

 si@surginno.co.uk

 www.surginno.com

 [@SurgicalInnovationsLtd](https://www.linkedin.com/company/SurgicalInnovationsLtd)

 [@SurgInno](https://twitter.com/SurgInno)

 [@SurgicalInnovations](https://www.youtube.com/SurgicalInnovations)

Contact Austria:

HCP Healthcarepartner GmbH
Wolfgangsee Straße 26, 5023 Salzburg-Guggenthal, Österreich
T +43 (0)662/62 12 43 • F -11 • office@hcp-austria.com

hcp-austria.com