

# ADVANCING A+ FOR LIFE

DISPOSABLE SKIN STAPLER

## OPEN INSTRUMENTS

TCYP Series



TONCARE  
Medical Instrument  
Manufacture Center



## ADVANCING FOR LIFE DISPOSABLE SKIN STAPLE

### Open Nose Piece and Low Profile

Open nose piece and low profile for clear view of staple placement. Arrow guide for accurate and consistent placement.

### Visible Staple Counter

Provide visible staple counter for constant indication of staple supply and for efficient use.

### Ergonomics Design

Ergonomics design of firing trigger for more stable firing.

### Staples Remover

Special staples remover can reduce the pressure on skin and relieve pain.



### Product Overview

TYCP Series has designed to deliver rectangular, stainless steel staples for routine wound closure.

Product Code	Closed Staple Height	Number Of Staples	Material Of Staples
TCYP(W)-35(A)	(W)7.0×(H)4.5	35	Medical stainless steel



- Streamlined handle design for comfortable and convenient operation.
- Small skin wound to avoid pain of traditional suture clearing.
- Smooth staple formation, accurate and safe.
- Save suture time and as far as possible to reduce pain from patient.
- Sterile medical products, single use, wide application scope.