

Experts for medical co-operations.

hcp+



THOR LLETZ LOOP ELECTRODES®

The new name for quality and precision
in electrosurgical instruments



THOR  **ELECTRODES®**

hcp-austria.com

LLETZ Loop Electrodes



0410 LLETZ SQUARE 10mm x 7mm (13cm)



0420 LLETZ LOOP 10mm x 5mm (13cm)



0430 LLETZ LOOP 15mm x 6mm (13cm)



0440 LLETZ LOOP 20mm x 10mm (13cm)



0450 LLETZ LOOP 10mm x 10mm (13cm)



FEATURES & BENEFITS



Nine different sizes available.
Colour-coded and disposable.



0460 LLETZ LOOP 15mm x 12mm (13cm)



0470 LLETZ LOOP 20mm x 12mm (13cm)



Accur
Generator



LLETZ LOOP 20mm x 15mm (13cm)

0480



LLETZ LOOP 20mm x 8mm (13cm)

0490



Fine Tungsten wire loop for better heat resistance. Full insulation throughout shaft.



Easy to use, strong 'Y' construction.



0009

5 inch (12.7cm) Ball Electrode



0035

Push Button Pencil with THOR PTFE Blade



0037

Push Button Pencil with Modified THOR PTFE Blade



+44 (0)7498 850552



admin@waves-medical.com



www.waves-medical.com



Waves Medical Office:

167-169 Great Portland Street
5th Floor
London
W1W 5PF

Additional Items



The Accur 400P is the most advanced ESU in our range. Delivering up to 400 Watts of power it can perform simultaneous coagulation on two active electrodes at the same time.

CONTACT AUSTRIA:

HCP Healthcarepartner GmbH

Wolfgangsee Straße 26, 5023 Salzburg-Guggenthal, Österreich
T +43 (0)662/62 12 43 • F -11 • office@hcp-austria.com



The Secur 200 and Secur 400 ESU's deliver power, performance and great functionality in a compact unit. They come with a 5.7" touchscreen for easy navigation and control.