

Experts for medical co-operations.

hcp+

Healthcare Products



THOR LAPAROSCOPIC ELECTRODES[®]

The new name for quality and precision
in electrosurgical instruments



THOR  **ELECTRODES[®]**

hcp-austria.com

PTFE Coated Laparoscopic Electrodes



- 0017 THOR PTFE Laparoscopic Blade 33cm
- 0017S THOR PTFE Laparoscopic Blade 33cm Split Stem



- 0018 THOR PTFE Laparoscopic Spatula 33cm
- 0018S THOR PTFE Laparoscopic Spatula 33cm Split Stem



- 0019 THOR PTFE Laparoscopic Curved Blade 33cm
- 0019S THOR PTFE Laparoscopic Curved Blade 33cm Split Stem



- 0020 THOR PTFE Laparoscopic L-Hook 33cm
- 0020S THOR PTFE Laparoscopic L-Hook 33cm Split Stem
- 0020X THOR PTFE Laparoscopic L-Hook 45cm
- 0020XS THOR PTFE Laparoscopic L-Hook 45cm Split Stem

FEATURES & BENEFITS



THOR Laparoscopic Electrodes are coated in high-grade PTFE, making surgery faster and less disrupted.



- 0021 THOR PTFE Laparoscopic J-Hook 33cm
- 0021S THOR PTFE Laparoscopic J-Hook 33cm Split Stem
- 0021X THOR PTFE Laparoscopic J-Hook 45cm
- 0021XS THOR PTFE Laparoscopic J-Hook 45cm Split Stem



- 0022 THOR PTFE Laparoscopic Needle 33cm
- 0022S THOR PTFE Laparoscopic Needle 33cm Split Stem





0024 THOR PTFE Laparoscopic Ball 33cm

0024S THOR PTFE Laparoscopic Ball 33cm Split Stem



0027 THOR PTFE Laparoscopic Thin Curved 33cm

0027S THOR PTFE Laparoscopic Thin Curved 33cm Split Stem



0100 THOR PTFE Laparoscopic Wire L-Hook 33cm

0100S THOR PTFE Laparoscopic Wire L-Hook 33cm Split Stem

0100X THOR PTFE Laparoscopic Wire L-Hook 45cm

0100XS THOR PTFE Laparoscopic Wire L-Hook 45cm Split Stem



0105 THOR PTFE Laparoscopic Wire J-Hook 33cm

0105S THOR PTFE Laparoscopic Wire J-Hook 33cm Split Stem



Reintroduction during surgery is minimised whilst lower power settings make for safer surgery.



Disposable electrodes provide a new fully insulated shaft for each case, preventing stray current and unintended burns.



Waves Medical Generator



Standard End - Compatible with standard pencils



Split Stem - Compatible with standard foot pedals



+44 (0)7498 850 552



admin@waves-medical.com



www.waves-medical.com



Waves Medical Office:

Waves Medical Ltd
167-169 Great Portland
Street, 5th Floor LONDON
W1W 5PF
UNITED KINGDOM

CONTACT AUSTRIA:

HCP Healthcarepartner GmbH

Wolfgangsee Straße 26, 5023 Salzburg-Guggenthal, Österreich

T +43 (0)662/62 12 43 • F -11 • office@hcp-austria.com