



SINGLE USE TRANSVERSE CUTTING LINEAR STAPLER AND RELOADS




SUZHOU LANTEX MEDICAL TECHNOLOGY COMPANY LIMITED

Building 7, 58 Xinting Road, High-tech Industrial Zone, 215151 Suzhou, China

Tel: +86 512 6616 7488

Fax: +86 512 6616 7488





PRODUCT SPECIFICATIONS

PRODUCT CODE	CARTRIDGE COLOR	STAPLE ROWS	STAPLE CROWN LENGTH (mm)	STAPLE LEG LENGTH (mm)	CLOSED STAPLE HEIGHT (mm)	STAPLING LINE LENGTH (mm)	STAPLE COUNT
TCLS40W	White	4	3.0	2.5	1.0	43	42
TCLS40B	Blue	4	3.0	3.8	1.5	43	42
TCLS40G	Green	4	3.0	4.5	2.0	43	42

LTX2017TCLS01EN